

4D Vaginal Hifu And Face Hifu System

Item Number: BL-509



Introduction

Professional 2-in-1 vaginal HIFU and 4D face HIFU system for professional medical aesthetic clinics. Non-invasive skin tightening, wrinkle removal, body contouring, and vaginal rejuvenation with 4D multi-line and vaginal HIFU probes. Stimulates collagen, no downtime, safe and effective.

[Learn More](#)

Application	Description	Key Benefit
Facial Wrinkle Reduction	Utilizes the DS-1.5 and DS-3.0 cartridges to deliver focused ultrasound to the epidermis and dermis, stimulating collagen and elastin production in areas such as crow's feet, forehead lines, and perioral wrinkles.	Smooths fine lines and improves skin texture with a natural, no-downtime treatment.
Full-Face Lifting	Deploys the DS-4.5 cartridge at 4.5 mm depth to target the SMAS layer, lifting sagging cheeks, jowls, and brows for a non-surgical facelift effect.	Provides a visible lift and tightened facial contour without incisions.
Jawline and Chin Contouring	Uses the 4.5 mm and optional 8 mm cartridges to tighten submental skin and redefine the jawline, reducing double chin appearance.	Creates a sharper, more youthful facial profile.
Body Contouring and Fat Reduction	Optional 13 mm cartridge focuses thermal energy into subcutaneous fat layers, inducing lipolysis and reducing localized fat deposits on the abdomen, flanks, and thighs.	Non-invasive body sculpting with measurable circumference reduction.
Acne Scar and Pore Refinement	Shallow 1.5 mm and 3.0 mm probes stimulate collagen remodeling to fill depressed scars and minimize enlarged pores.	Restores smooth, even skin surface and reduces pore size.
Skin Tone Improvement & Brightening	The 3.0 mm cartridge accelerates cell turnover and collagen synthesis, leading to improved skin tone, reduced pigmentation, and a more radiant complexion.	Achieves youthful luminosity and even skin tone.
Vaginal Tightening	Vaginal-specific 3.0 mm and 4.5 mm cartridges deliver HIFU to the lamina propria and muscular layer, causing collagen denaturation and subsequent tightening of the vaginal canal.	Enhances vaginal tone, improves intimate comfort, and supports sexual wellness.
Vaginal Lubrication and Sensitivity Enhancement	Submucosal 3.0 mm treatment accelerates cellular regeneration and promotes natural lubrication by revitalizing glandular function.	Alleviates vaginal dryness and heightens sensitivity for improved intimate satisfaction.
Eye Bag and Dark Circle Reduction	Precision targeting of the delicate periorbital area with the 1.5 mm cartridge firms skin and reduces the appearance of under-eye bags and dark circles.	Refreshes the eye area, making it look brighter and more youthful.

Parameter	Specification
Model	BL-509
Modes	Vaginal HIFU and Facial 4D HIFU
Power	400W
Voltage	100-240V, 50Hz/60Hz
Energy	Up to 2.0 J
Vaginal Cartridges	3.0 mm and 4.5 mm
Facial Cartridges	1.5 mm, 3.0 mm, 4.5 mm (8 mm and 13 mm optional)
Shots per Cartridge	10,000 shots
Total Shot Capacity	50,000 shots

Parameter	Specification
Facial Probe DS-4.5	4 MHz, penetration depth 4.5 mm (targets SMAS layer)
Facial Probe DS-3.0	8 MHz, penetration depth 3.0 mm (targets dermal collagen layer)
Facial Probe DS-1.5	10 MHz, penetration depth 1.5 mm (targets epidermal tissue)
Vaginal Probe 4.5 mm	Delivers energy to SMAS layer for contraction and tightening
Vaginal Probe 3.0 mm	Delivers energy to submucosal tissue for collagen regeneration
Gross Weight	25 kg
Package Dimensions	65 x 40 x 70 cm